

Bars & Frameworks

customer manual

Table of contents

| | |
|---------------------------------|---|
| Bars & frameworks overview | 1 |
| Bars & frameworks designs | 2 |
| Bars & frameworks design matrix | 3 |
| Impression components | 4 |

Bars & frameworks overview

Bars & frameworks introduction

CAD/CAM Patient Specific Restorations are the future of restorative implant dentistry and BioHorizons is leading the way with Bars and Frameworks. These superstructures provide simplified laboratory procedures for implant overdentures, fixed-hybrid prostheses and fixed prostheses.

With precision difficult to match using conventional laboratory techniques, BioHorizons Bars and Frameworks are one-piece milled titanium alloy structures with a passive fit. One-piece milled frameworks are significantly stronger than cast bars and do not require soldering or welding. The result is a durable restoration with a precise fit.


BioHorizons validated milling centers design bars and frameworks in CAD from a laboratory or clinician submitted order form and then send the design electronically to the customer for design verification. After verification, the CAD design is milled into a bar or framework. Then it's ready for overdenture processing or the addition of porcelain or acrylic resin with minimal finishing required by the customer.

Copy-milled design

Laboratory technicians can create their own unique bar or framework design with a resin pattern and send it to a BioHorizons validated milling center with the master cast to be scanned and milled. Using a copy-milling technique, the milling center creates a one-piece, titanium alloy replica of the design provided, to fit the master cast.

BioHorizons bars & frameworks

- Designed passive fit
- No soldered or welded joints
- Lightweight
- CAD/CAM precision
- No capital investment
- Laboratory design control
- No waxing and casting
- Biocompatible alloy




Bars & frameworks designs

Design options

BioHorizons Bars and Frameworks are compatible with all internal implant systems and match a wide variety of implant and abutment connections.


Ellipse




Dolder




C+M Dolder Macro



C+M Eggshape Micro




C+M Eggshape Macro




Primary Bar Straight




Primary Bar 2° Tapered




Primary Bar 4° Tapered




Round 3.0mm



Hader Bar 3mm



Fixed Implant Denture



Freeform

As specified by the prescribing clinician within controlled design limitations

Bars & frameworks design matrix

Implant level bars


| | Description | Minimum | Maximum |
|---|----------------------------------|---------|----------------|
| A | Platform Seating Diameter | 3.0mm | 5.7mm |
| B | Total Cylinders | 2 | 10 |
| C | Bar Span Between Cylinders | 0mm | 30mm |
| D | Bar Height incl. Cylinder | 2.5mm | 11mm |
| E | Maximum Angulation Per Cylinder* | 0° | 20° |
| F | Cylinder Diameter | 3.0mm | 10.0mm |
| G | Distal Extension | 0mm | 1.5x AP Spread |

*Does not apply to third-party compatible implant lines.

Abutment level bars

| | Description | Minimum | Maximum |
|---|---------------------------------|---------|----------------|
| A | Platform Seating Diameter | 4.8mm | 7mm |
| B | Total Cylinders | 2 | 10 |
| C | Bar Span Between Cylinders | 0mm | 30mm |
| D | Bar Height incl. Cylinder | 2.5mm | 11mm |
| E | Maximum Angulation Per Cylinder | 0° | 20° |
| F | Cylinder Diameter | 4.8mm | 10.0mm |
| G | Distal Extension | 0mm | 1.5x AP Spread |
| H | Screw Channel Angulation | 0° | 25° |

Note: BioHorizons reserves the right to reject any case request that does not meet regulatory requirements.



Impression components



Multi-unit CAD/CAM Bar Screws

| | |
|------------------|--|
| PXMUPSA | Multi-unit CAD/CAM Bar Screw |
| PXMUPSA25 | Multi-unit CAD/CAM Bare SCrew (pack of 25) |

For attaching CAD/CAM bar restorations to the Multi-unit abutments. Titanium alloy. Hand-tighten or torque to 15Ncm with .050" (1.25mm) Hex Driver or Ballpoint Driver (CBDRL), depending on application. The Ballpoint Driver is required for bars with an angled screw channel.

The Multi-unit CAD/CAM bar screw will fixate a CAD/CAM Bar with either straight or angled screw channels to multi-unit abutments for currently marketed abutment level CAD/CAM bars.



Abutment Screws

| | |
|---------------|-----------------------------|
| PXAS | Abutment Screw |
| PXAS25 | Abutment Screw (pack of 25) |

PXAS fixates internal hex, implant level bars. Final torque: 30 Ncm with .050 (1.25mm) Hex Driver.

Direct Offices

BioHorizons USA
888-246-8338 or 205-967-7880

BioHorizons Canada
866-468-8338

BioHorizons Spain
+34 91 713 10 84

BioHorizons UK
+44 (0)1344 752560

BioHorizons Chile
+56 (2) 23619519

BioHorizons Italy
800-063-040

BioHorizons Mexico
800-953-0498

Distributors

For contact information in our 90 countries, visit biohorizons.com



shop online at
store.biohorizons.com

BioHorizons®, Laser-Lok®, MinerOss®, AutoTac®, Mem-Lok®, TeethXpress®, IntraSpin®, L-PRF® and Xpression® are registered trademarks of BioHorizons Implant Systems Inc. Genate™ is a trademark of BioHorizons Implant Systems Inc. Unigrip™ is a trademark of Nobel Biocare AB. Zimmer® Dental ScrewVent® and Tapered ScrewVent® are registered trademarks of Zimmer, Inc. AlloDerm SELECT™, AlloDerm SELECT GBR™ and NovoMatrix™ are trademarks of Allergan, an AbbVie company. Grafton® DBM is a registered trademark of Medtronic, Inc. Cytoplast® is a registered trademark of Osteogenics Biomedical, Inc. Puros Dermis is a registered trademark of Zimmer Biomet. Mucograft is a registered trademark of Ed. Geistlich Sogne Ag Fur Chemische Industrie. Symbios PerioDerm is a registered trademark of Dentsply Sirona. Hu-Friedy® is a registered trademark of Hu-Friedy Mfg. Co., LLC. Spiralock® is a registered trademark of Spiralock Corporation. Pomalux® is a registered trademark of Westlake Plastics Co. Locator® is a registered trademark of Zest Anchors, Inc. Delrin® is a registered trademark of E.I. du Pont de Nemours and Company. Bio-Gide® is a registered trademark of Edward Geistlich Sohne AG Fur Chemische Industrie. BioMend® is a registered trademark of Zimmer Biomet Dental. Not all products shown or described in this literature are available in all countries. We are proud to be registered to ISO 13485:2016, the international quality management system standard for medical devices, which supports and maintains our product licences with Health Canada and in other markets around the globe. Original language is English. ©BioHorizons. All Rights Reserved.

L02042
L02042



REV B MAR 2025
REV B MAR 2025