

Elos Accurate for Dentsply Sirona 3Shape Library Rev D Release Notes

What's New

- Renamed library name from “Astra Tech” to “Dentsply Sirona”
- New abutment kits: Ankylos Balance Base and Astra Uni Abutment.
- Added non-engaging hybrids (HBN) to Ankylos Balance Base and Astra Uni Abutment
- Initial release of tall non-engaging hybrids (HBN-H) for Ankylos Balance Base and Astra Uni Abutment

Supported Components



Abutments

- Custom Ti for Astra TX, Astra EV
- Hybrid Base for Astra TX (HBE, HBE-H)
- Hybrid Base for Astra EV (HBE, HBE-H)
- Hybrid Base for Ankylos Balance Base (HBN, HBN-H)
- Hybrid Base for Astra TX UNI Abutment (HBN, HBN-H)
- Hybrid Base for Astra EV UNI Abutment (HBN, HBN-H)

Scan Bodies

- Elos I-O (TX, EV, TX UNI Abutment, EV UNI Abutment)
- 3Shape (TX, EV)
- NT Trading (TX)

Analogs

- Elos Analogs (PMA)

Workflows

- 28° Angled Screw Channel with Hybrid Abutments
- Custom Ti, Temporary & Zirconia restorations

