

Library Name: V-DME010 – Elos for Nobel

Revision: A

Restorations: *(Check all that apply)*

Custom Titanium Hybrid Restoration Printed Model Bars

Changes made: *(list all changes made to library)*

1. Includes Custom Ti for Nobel Active and Nobel Replace previously within V-DME008 Rev D.
2. Added all files for Nobel Active 3.0, NP, RP, and WP Hybrid Abutments. The 3.0 connection just has an engaging hybrid and the NP, RP, WP all have engaging and non-engaging hybrid abutments. Scan body options include both Elos scan bodies and 3Shape scan bodies
3. Added all files for NobelReplace NP, RP, WP, and 6.0 Hybrid Abutments (engaging and non-engaging). This includes scan body options for Elos scan bodies and 3Shape scan bodies
4. Added all files for Nobel Multi-Unit Non-Engaging Hybrid Abutments. Elos scan bodies only.
5. Added VCD Elos Hybrid Nobel Active, VCD Elos Hybrid Nobel Replace and VCD Elos Hybrid Nobel MU to the Vulcan – Elos Hybrid Base system category
6. Added all files for Branemark NP, RP, and WP Custom Titanium abutments. This includes scan bodies for Elos and 3Shape.
7. Added VCD Elos Ti Branemark to the Vulcan – Elos Custom Ti system category.