

1 Identification

- **Product identifier**
- **Trade name:** **NexxZr T/ + / S (Disc)**
- **Application of the substance / the mixture**
Manufacture of dental prothesis
Yttrium-stabilized zirconium oxide disc for the CAD/CAM technology
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Manufacturer:
Sagemax Bioceramics, Inc.
34210 9th Avenue South, Suite 118
Federal Way, WA 98003 USA
Tel. +1 253-214-0389 / +1 877-386-0389
Fax +1 888-856-2615

- EC Representative:
AB Ardent
Generatorgatan 8
SE-19560 Arlandastad / Sweden
Tel. +46 8594 412 57
Fax +46 8591 161 83

- **Information department:**
Quality Assurance / Regulatory Affairs
USA: support@sagemax.com
EU: info@sagemax.com
- **Emergency telephone number:**
24 Hour Emergency Assistance:
Emergency-Call USA - Infotrac: 1-800-535-5053
Emergency-Call Canada - Canutec: 1-613-996-6666

- General SDS Assistance:
US: 1-800-533-6825
Canada: 1-800-263-8182

2 Hazard identification

- **Classification of the substance or mixture**
The product is not classified, according to the Globally Harmonized System (GHS).
-
- **Label elements**
 - **GHS label elements** Void
 - **Hazard pictograms** Void
 - **Signal word** Void
 - **Hazard statements** Void
 - **Classification system:**
 - **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**
- | | | |
|------------|---|----------------|
| HEALTH | 0 | Health = 0 |
| FIRE | 0 | Fire = 0 |
| REACTIVITY | 0 | Reactivity = 0 |

Safety Data Sheet

according to HPR, Schedule 1

Printing date 04/04/2019

Version number 1

Reviewed on 04/04/2019

Trade name: **NexxZr T/ + / S (Disc)**

(Contd. of page 1)

· **Other hazards** Avoid breathing sanding dust /swarf.

3 Composition/Information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** Mixture of the substances listed below with nonhazardous additions.

- **Dangerous components:**

CAS: 1314-23-4	zirconium oxide	60-100% w/w *
CAS: 1314-36-9	yttrium oxide	1-5% w/w *
CAS: 12055-23-1	hafnium dioxide	1-5% w/w *

* Actual concentration ranges are withheld as a trade secret.

4 First aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:** Grinding dust: Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:**
Grinding dust: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:**
Grinding dust: Rinse out mouth and then drink plenty of water.
If symptoms persist consult doctor.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Firefighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:** The product is not flammable.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** No special measures required.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

CA

(Contd. on page 3)

Safety Data Sheet

according to HPR, Schedule 1

Printing date 04/04/2019

Version number 1

Reviewed on 04/04/2019

Trade name: **NexxZr T/ + / S (Disc)**

(Contd. of page 2)

7 Handling and storage

- **Handling:**
- **Precautions for safe handling**
Only adequately trained personnel should handle this product.
For use in dentistry only.
Ensure good ventilation/exhaustion at the workplace.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store only in the original receptacle.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s)** No further relevant information available.

8 Exposure controls/ Personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**

- **Components with limit values that require monitoring at the workplace:**

CAS: 1314-23-4 zirconium oxide

EL	Short-term value: 10 mg/m ³ Long-term value: 5 mg/m ³ as Zr
----	---

CAS: 1314-36-9 yttrium oxide

EL	Long-term value: 1 mg/m ³ as Y
----	--

CAS: 12055-23-1 hafnium dioxide

EL	Long-term value: 0.5 mg/m ³ as Hf
----	---

- **Additional information:** The lists that were valid during the creation were used as basis.

- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**
Usual hygienic measures for dental practice and dental laboratories.
Keep away from foodstuffs, beverages and feed.
Wash hands before breaks and at the end of work.
Avoid breathing sanding dust /swarf.
- **Breathing equipment:**
Mechanical processing:
Use respiratory protective device against the effects of dust.
- **Recommended filter device for short term use:** Filter P1
- **Protection of hands:**
Protective gloves should always be worn during mechanical and thermal processing.
After use of gloves apply skin-cleaning agents and skin cosmetics.
- **Material of gloves**
Mechanical processing:
Leather gloves
Strong gloves
- **Penetration time of glove material -**

(Contd. on page 4)

CA

Safety Data Sheet

according to HPR, Schedule 1

Printing date 04/04/2019

Version number 1

Reviewed on 04/04/2019

Trade name: **NexxZr T/ + / S (Disc)**

(Contd. of page 3)

- **Eye protection:**
Always wear safety goggles during mechanical processing (grinding, sawing /cutting, drilling, milling).

9 Physical and chemical properties

- **Information on basic physical and chemical properties**

- **General Information**

- **Appearance:**

· Form:	Solid
· Color:	Whitish
· Odor:	Odorless
· Odor threshold:	Not determined.

- **pH-value:** Not applicable.

- **Change in condition**

· Melting point/Melting range:	Undetermined.
· Boiling point/Boiling range:	Undetermined.

- **Flash point:** Not applicable.

- **Flammability (solid, gaseous):** Product is not flammable.

- **Auto igniting:** Product is not selfigniting.

- **Danger of explosion:** Product does not present an explosion hazard.

- **Vapor pressure:** Not applicable.

· Density at 20 °C:	3-3.5 g/cm ³
· Relative density	Not determined.
· Vapor density	Not applicable.
· Evaporation rate	Not applicable.

- **Solubility in / Miscibility with**

· Water:	Insoluble.
-----------------	------------

- **Partition coefficient (n-octanol/water):** Not determined.

- **Viscosity:**

· Dynamic:	Not applicable.
· Kinematic:	Not applicable.

- **Other information** No further relevant information available.

10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability** Stable under normal handling and storage conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** None under normal conditions of storage and use.

CA

(Contd. on page 5)

Safety Data Sheet

according to HPR, Schedule 1

Printing date 04/04/2019

Version number 1

Reviewed on 04/04/2019

Trade name: **NexxZr T/ + / S (Disc)**

(Contd. of page 4)

11 Toxicological information

- **Information on toxicological effects**
- **Additional toxicological information:**
- **Carcinogenic categories**
- **NTP (National Toxicology Program)**

None of the ingredients is listed.

12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:** Not hazardous for water.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**
Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

14 Transport information

- | | |
|---------------------------------------|-----------------|
| · UN-Number | Void |
| · DOT, TDG, ADN, IMDG, IATA | Void |
| · UN proper shipping name | Void |
| · DOT, TDG, ADN, IMDG, IATA | Void |
| · Transport hazard class(es) | Void |
| · DOT, TDG, ADN, IMDG, IATA | Void |
| · Class | Void |
| · Packing group | Void |
| · DOT, TDG, IMDG, IATA | Void |
| · Environmental hazards: | Not applicable. |
| · Special precautions for user | Not applicable. |

(Contd. on page 6)

CA

Safety Data Sheet

according to HPR, Schedule 1

Printing date 04/04/2019

Version number 1

Reviewed on 04/04/2019

Trade name: NexxZr T/ + / S (Disc)

(Contd. of page 5)

- | | |
|--|--|
| · Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code | Not applicable. |
| · Transport/Additional information: | Product is not classified as a dangerous good for transport (ADR, IMDG, IATA). |
| · UN "Model Regulation": | Void |

15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Sara**
- **Section 355 (extremely hazardous substances):**
None of the ingredients is listed.
- **Section 313 (Specific toxic chemical listings):**
None of the ingredients is listed.
- **TSCA (Toxic Substances Control Act):**
All ingredients are listed.
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **National regulations:**
- **Other regulations, limitations and prohibitive regulations**
The product is a medical device according to the Directive 93/42/EEC.
This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.
- **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Date of the latest revision of the safety data sheet** 04/04/2019 / -
- **Abbreviations and acronyms:**
 IMDG: International Maritime Code for Dangerous Goods
 DOT: US Department of Transportation
 IATA: International Air Transport Association
 EINECS: European Inventory of Existing Commercial Chemical Substances
 ELINCS: European List of Notified Chemical Substances
 CAS: Chemical Abstracts Service (division of the American Chemical Society)
 NFPA: National Fire Protection Association (USA)
 HMIS: Hazardous Materials Identification System (USA)
 PBT: Persistent, Bioaccumulative and Toxic
 vPvB: very Persistent and very Bioaccumulative

CA