

## BioHorizons Trimmable Prosthetics 3Shape Library

### Rev 03 Release Notes

### What's New

- Added Multi-unit Contour Cover Cap (PXMUCCC) to BioHorizons Trimmable Prosthetics library.

### Supported Components



#### Abutments

- Easy Ti Temporary Abutment (Hexed, Non-Hexed)
- PEEK Temporary Cylinder (Hexed, Non-Hexed)
- BioHorizons Multi-Unit Ti Coping

#### Scan Bodies

- BioHorizons: PEEK Scan Abut, Snap Scan (3 heights), Multi-unit scan caps\* (e.g., PXMUCCC)
- 3rd party: Elos\*, 3Shape, NT Trading, Anatotemp SC

#### Analogs

- BioHorizons authentic lab analogs (e.g., PYIA)

#### Workflows

- Full-Arch Restorations
- Denture Conversion

**\*NOTE:** BioHorizons Multi-Unit Ti Coping only supports indicated scan bodies.

