

# BioHorizons Trimmable Prosthetics 3Shape Library

## Rev 04 Release Notes

### What's New

- Added new titanium scan bodies (TP3TSB, PYTSB, PGTSB, PBTBS and PXMUTSB)
- Added direct pick-up kits back into the library

### Supported Components



#### Abutments

- Easy Ti Temporary Abutment (Hexed, Non-Hexed)
- PEEK Temporary Cylinder (Hexed, Non-Hexed)
- BioHorizons Multi-Unit Ti Coping

#### Scan Bodies

- BioHorizons: Titanium Scan bodies\*, PEEK Scan Abut, Snap Scan (3 heights), Multi-unit scan caps\* (e.g., PXMUCCC)
- 3rd party: Elos\*, 3Shape, NT Trading

#### Analog

- BioHorizons authentic lab analogs (e.g., PYIA)

#### Workflows

- Full-Arch Restorations
- Denture Conversion

**\*NOTE:** BioHorizons Multi-Unit Ti Coping only supports indicated scan bodies.

