

## BioHorizons CAD-CAM & BioHorizons Analogs ExoCAD Library Rev F Release Notes

### What's New

- Added Multi-unit Contour Cover Cap (PXMUCCC) to BioHorizons Hybrid Non-Hexed library.

### Supported Components

exocad

#### Abutments

- Hybrid Base (Hexed, Non-Hexed, 2mm, Tissue Level)
- Laser-Lok Ti Base
- Custom Ti Abutment (with and without Laser-Lok)

#### Scan Bodies

- BioHorizons: PEEK Scan Abut, Snap Scan (3 heights), Multi-unit scan caps (e.g., PXMUCCC)
- 3<sup>rd</sup> party: Elos, 3Shape, NT Trading, Anatotemp SC

#### Analogs

- Elos Accurate analogs – BioHorizons compatible (e.g., PMA-BHI35-1)
- BioHorizons authentic lab analogs (e.g., PYIA)

#### Workflows

- 15° Angled Screw Channel with Hybrid Base Abutments (excluding 2mm and Tissue Level)
- Custom Ti, Temporary & Zirconia restorations

