

## Elos Accurate 3Shape Library (Rev D) Straumann - Neodent

### What's New

- Now supporting Neodent with Custom Ti

### Supported Components



#### Abutments

- Custom Ti for Straumann Bone Level, Neodent
- Hybrid Base for Straumann Bone Level (HBE, HBE-H, HBN)
- Hybrid Base for Straumann Screw-retained (HBN, HBN-H)

#### Scan Bodies

- Elos I-O (Bone Level, Screw-retained, Neodent)
- 3Shape (Bone Level)
- NT Trading (Bone Level)
- Anatotemp SC (Bone Level RC)

#### Analogs

- Elos Analogs (PMA)

#### Workflows

- 28° Angled Screw Channel with Hybrid Abutments
- Custom Ti, Temporary & Zirconia restorations

